## Quiz 5

## TOTAL POINTS 16

1.	Which shows the highest signal intensity in T1 -weighted image?	1 point
	T1 = 200 ms, T2 = 100 ms	
	T1 = 1000 ms, T2 = 500 ms	
	T1 = 1000 ms, T2 = 100 ms	
	T1 = 300 ms, T2 = 50 ms	
2.	Tissue with short T1 shows ( A ) signal intensity on T1 -weighted image and tissue with short T2 shows ( B ) signal intensity on T2 -weighted image. (A , B)	1 point
	high, low	
	O high, high	
	O low, low	
	O low, high	
3.	When FOV increases, the maximum k-space range also increases.	1 point
	○ True	
	False	
4.	The step size in k-space is inversely related to resolution in image.	1 point
	○ True	
	<ul><li>False</li></ul>	
5.	Which is an essential element to correct magnetic field inhomogeneity in spin echo imaging?	1 point
	180° RF pulse	
	Phase encoding gradient	
	○ 90° RF pulse	
	Slice select gradient	
6.	Which of the follows is an intrinsic tissue-specific factor to determine contrast?	1 point
	Proton density	
	O TE	
	○ Flip angle	
	○ TR	

7.	Wh	ich statement is correct?		1 point
	0	Scan time is not related to TR.		
	0	Tissue with high proton density shows low signal intensity in proton density-weighted images.		
	0	To reduce T1 contrast, TR should be minimized.		
	•	T2 -weighted images are sensitive to TE.		
		8. In gradient echo imaging, what is optimal flip angle (Ernst angle) for a tissue with T1=1sec when TR = 100 msec?	1 point	
		○ 15°		
		○ 45°		
		○ 35°		
		Select the best TE value that maximizes T2 contrast between two tissues with T2 values of 60 msec and 100 msec.	1 point	
		<ul><li>● 80 msec</li></ul>		
		30 msec		
		As short as possible		
		150 msec		
		10. Select the best TR value that maximizes proton density contrast between two tissues with T1 values of 800 msec and 1000 msec.	1 point	
		● 600 msec		
		900 msec		
		1200 msec		
		○ 3000 msec		
		11. Select the best TE value that maximizes T1 contrast between two tissues with T1 values of 800 msec and 1000 msec and T2 values of 60 msec and 100 msec.	1 point	
		As short as possible		
		○ 80 msec		
		900 msec		
		○ 3000 msec		
		12. Brain lesion that contains magnetically various components shows low signal intensity in gradient echo imaging.	1 point	
		False		
		13. Which statement is true?	1 point	
		When sampling rate and number of points are constant, high frequency encoding gradient makes resolution poor.		
		When FOV increases, the frequency encoding gradient strength also increases for a fixed sampling rate.		
		The larger FOV, the better resolution		

 $\begin{tabular}{ll} \hline & When FOV increases, sampling rate also increases for a fixed frequency encoding gradient. \\ \hline \end{tabular}$ 

14.	The readout prephase gradient is used to generate gradient echo in gradient echo imaging, therefore it is not required spin echo imaging.	in 1	point
	○ True		
	<ul><li>False</li></ul>		
15	5. K-space region is sampled at interval of $\Delta$ Kx within the range of [-Kx $_{max}$ / 2, Kx $_{max}$ / 2] along the horizontal direction. A new k-space region is sampled at interval of $\Delta$ Kx $_{max}$ / 2 within range of [-Kx $_{max}$ / 2, Kx $_{max}$ / 2]. How are the fied of view (FOV) and the resolution affected in the new image compared to the first image?	1 point	
	Increasing FOV, decreasing resolution.		
	No change in FOV, increasing resolution		
	O Decreasing FOV, increasing resolution.		
	Increasing FOV, no change in resolution		
16	6. Which is the appropriate way of obtaining high resolution images?	1 point	
	High frequency encoding gradient (assuming that sampling rate and number of points remain constant)		
	O Increasing sampling rate (assuming that frequency encoding gradient and number of points remain constant)		
	○ Increasing FOV		
	O Decreasing matrix size (i.e., the number of the points)		

Answer of Q.10: 1200 msec

Answer of Q.13 When FOV increases, sampling rate also increases for a fixed encoding gradient.